



Configuration #GP2AP-SHF510-40EG
Elevating Table with Floor Mount Rotate Tube Stand
with 50 kW Generator and 40EG Dual Digital Panels

SHF 510AP-GI	50 kW, 500 mA, 500 mAs High Frequency Generator <ul style="list-style-type: none">• 500 mA., 125 kVp (500 mA. @100 kVp)• Eighteen mA stations, 10, 12.5, 16, 20, 25, 32, 40, 50, 64, 80, 100, 125, 160, 200, 250, 320, 400, 500• 0.001 to 10 second timing• kVp selection, 40 to 125 kVp in 1 kV increments• Digital displays• Automatic line compensation• APR, Anatomical Programming• AEC Electronics and (2) ion chambers• Generator integration with acquisition monitor• Exposure hand switch• UL/CE/FDA Listed• For Operation on 1 Phase 220/240VAC, see option for 3 phase
S-HL	Elevating, Four Way Float Top Table -750 Lb. Patient Load <ul style="list-style-type: none">• 86.5" x 31.5" float table top• 43.3" longitudinal and 9.5" lateral table top travel• Collision protection
RECEPTOR/ ACQUISITION GRIDS	Canon CXDI-40EG Fixed DR Panels, (Qty. 2) <ul style="list-style-type: none">• 17" x 17" (43cm x 43cm) imaging area• LANMIT 6, Amorphous Silicon and Scintillator• Gadolinium Oxysulfide Scintillator• 2688 x 2688 pixel array (7.2 million pixels)• 12 bit/4096 gray scale contrast• Preview image in about 5 seconds• Power Supply and Ethernet interface• Dual mode landscape/portrait DR cabinets (Qty. 2)• 10:1 ratio grid for DR cabinets (Qty. 2)• DICOM 3.0 complaint• DICOM Modality Worklist• 19" Touch screen monitor with generator integration• Computer, 80 GB• Canon LANMIX MLT processing software• Keyboard• Mouse
S-FMT	Floor Mount Tube Stand <ul style="list-style-type: none">• 11.8" Transverse travel, 180° rotation, vertical travel 15.7" to 74.75"• Platform tube mount• Electromagnetic brakes for longitudinal, transverse, vertical & rotation• 86.6" floor track provides 58.3" longitudinal tube travel
MIL-WH	Wall Stand <ul style="list-style-type: none">• Floor to Wall Mount with magnetic locks (84" column)
E7884FX	X-Ray Tube- 300,000 Heat Unit, 90° cable arms, 150 kVp, 0.6-1.2mm focal spot, 3" anode, 14" target, LS
70-63000-6	Manual Collimator, Includes <ul style="list-style-type: none">• Molded plastic housing• Bucky centering light• Tape measurer
10526-25	Cables, Pair - 25' long



Configuration #GP2AP-SHF510-40EG Elevating Table with Floor Mount Rotate Tube Stand with 50 kW Generator and 40EG Dual Digital Panels

Generator Options:

PED	Pedestal for control console, replaces table mount
CP	Wall mount for control console, replaces table mount
SHF 520AP	50 kW, 125 kVp 3 phase, 220VAC
SHF 525AP	50 kW, 150 kVp 3 phase, 220VAC
SHF 545AP	50 kW, 150 kVp 3 phase, 480VAC
SHF 645AP	64 kW, 150 kVp 3 phase, 480VAC, replaces 50 kW
SHF 845AP	80 kW, 150 kVp 3 phase, 480VAC, replaces 50 kW
HSS-1T	High speed starter, 1 tube

Table Options:

AEC	(See Generator Option)
-----	------------------------

Acquisition/Receptor Options: See page 47 for additional upgrade options

RADAX	Rotating DR cabinet, replaces landscape/portrait cabinet
-------	----------------------------------------------------------

Wall Stand Options:

TE-11	Tilt mechanism
-------	----------------

Tube Options:

E7254FX	400,000 H.U., 150 kVp, 0.6-1.2mm tube, replaces E7884FX
---------	---------------------------------------------------------

E7869FX	600,000 H.U., 150 kVp, 0.6-1.2mm tube, replaces E7884FX
---------	---------------------------------------------------------

Grid Options:

1011783444	10:1 ratio, 178 line grid for table, replaces 103 line
------------	--------------------------------------------------------

1011784072	10:1 ratio, 178 line grid for wall stand, replaces 103 line
------------	-------------------------------------------------------------

Collimator Options:

SM-302A-2	Automatic collimator replaces, manual collimator
-----------	--------------------------------------------------

Assistance:

(Does not include travel and expenses)

AST-DR	Onsite assistance for DR panel
--------	--------------------------------

AST-APDR	Onsite applications training
----------	------------------------------

